→ BPL_TEST2_PID_Fedbatch_reg6 script with FMPy 0.3.15

The key library FMPy ver 0.3.15 is installed.

After the installation a small application BPL_TEST2_PID_Fedbatch_reg6 is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google
    No LSB modules are available.
    Distributor ID: Ubuntu
    Description: Ubuntu 22.04.2 LTS
    Release:
                  22.04
    Codename:
                   jammy
%env PYTHONPATH=
    env: PYTHONPATH=
!wget https://repo.anaconda.com/miniconda/Miniconda3-py310_23.1.0-1-Linux-x86_64.sh
!chmod +x Miniconda3-py310_23.1.0-1-Linux-x86_64.sh
!bash ./Miniconda3-py310_23.1.0-1-Linux-x86_64.sh -b -f -p /usr/local
sys.path.append('/usr/local/lib/python3.10/site-packages/')
    --2023-09-13 18:32:00-- https://repo.anaconda.com/miniconda/Miniconda3-py310 23.1.0-1-Linux-x86 64.s
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.131.3, 104.16.130.3, 2606:4700::6810:8203,
    Connecting to repo.anaconda.com (repo.anaconda.com) | 104.16.131.3 | :443... connected.
    HTTP request sent, awaiting response... 200 OK
    Length: 74403966 (71M) [application/x-sh]
    Saving to: 'Miniconda3-py310_23.1.0-1-Linux-x86_64.sh'
    242MB/s
                                                                    in 0.3s
    2023-09-13 18:32:00 (242 MB/s) - 'Miniconda3-py310 23.1.0-1-Linux-x86 64.sh' saved [74403966/74403966
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Downloading and Extracting Packages
    Downloading and Extracting Packages
    Preparing transaction: done
    Executing transaction: done
    installation finished.
!conda update -n base -c defaults conda --yes
```

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done
```

!conda --version !python --version

> conda 23.7.4 Python 3.10.13

!conda install -c conda-forge fmpy --yes # Install the key package

Preparing transaction: done Verifying transaction: done Executing transaction: done

!conda install matplotlib --yes

```
Preparing transaction: done
    Verifying transaction: done
    Executing transaction: done
#!conda install scipy --yes
#!conda install xlrd --yes
#!conda install openpyxl --yes
```

→ BPL_TEST2_Fedbatch setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:

- FMU BPL_TEST2_Fedbatch_linux_om_me.fmu
- Setup-file BPL_TEST2_Fedbatch_fmpy_explore.py

```
%%bash
git clone https://github.com/janpeter19/CONF_2023_10_MODELICA15
    Cloning into 'CONF_2023_10_MODELICA15'...
%cd CONF_2023_10_MODELICA15
    /content/CONF_2023_10_MODELICA15
run -i BPL_TEST2_PID_Fedbatch_reg6_fmpy_explore.py
    Linux - run FMU pre-comiled OpenModelica 1.21.0
    Model for bioreactor has been setup. Key commands:
     - par() - change of parameters and initial values
                   - change initial values only
     - init()


    simulate and plot

     - simu()
     - newplot() - make a new plot
                 - show plot from previous simulation
                   - display parameters and initial values from the last simulation
     - disp()
     - describe() - describe culture, broth, parameters, variables with values/units
    Note that both disp() and describe() takes values from the last simulation
    and the command process_diagram() brings up the main configuration
    Brief information about a command by help(), eg help(simu)
    Key system information is listed with the command system_info()
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
```

```
import warnings
warnings.filterwarnings("ignore")
```

▼ BPL_TEST2_Fedbatch - demo

process_diagram()



describe('culture'); print(); #describe('liquidphase')

Pump schedule pa

Simplified text book mode - only substrate S and cell concentration X

```
par(Y=0.40, qSmax=1.0, Ks=0.1)
init(V_0=1e3, VX_0=1e3, VS_0=10*1e3)
\# Culture parameters
# Process initialization
par(S_in=600)
par(t_start=4.3, F_start=4, mu_feed=0.2, F_max=35)
par(S_ref=0.1) # Substrate controller
par(t_regStart=4.3)
par(uMax=50)
newplot()
ax2.set_ylim([0, 0.45]); ax3.set_ylim([0, 1])
setLines(['-']);
par(K=30, Ti=0.5)
simu(8)
```

```
Fedbatch cultivation with substrate control
         10
       X [g/L]
          5
        0.4
       mu [1/h]
         0.2
        0.0
        1.0
      S [g/L]
S
        0.0
        <u>F</u>
disp(mode='long')
    bioreactor.V_0 : V_0 : 1000.0
    bioreactor.m 0[1] : VX 0 : 1000.0
    bioreactor.m_0[2] : VS_0 : 10000.0
    bioreactor.culture.Y : Y : 0.4
    bioreactor.culture.qSmax : qSmax : 1.0
    bioreactor.culture.Ks : Ks : 0.1
     feedtank.V_0 : feedtank_V_0 : 10.0
     feedtank.c_in[2] : S_in : 600
    dosagescheme.mu_feed : mu_feed : 0.2
    dosagescheme.F_0 : F_0 : 0.0
    dosagescheme.t_start : t_start : 4.3
    dosagescheme.F_start: F_start: 4
    dosagescheme.F_max : F_max : 35
     substrateSensor.x_0 : Sensor_x_0 : 0
    substrateRef.k : S_ref : 0.1
    t_regStart : t_regStart : 4.3
    PIDreg.K : K : 30
    PIDreg.Ti : Ti : 0.5
    PIDreg.I 0 : I 0 : 0
    uMax : uMax : 50
describe('mu')
    Cell specific growth rate variable : 0.199 [ 1/h ]
describe('parts')
     ['bioreactor', 'bioreactor.culture', 'dosagescheme', 'feedtank', 'PIDreg', 'substrateRef', 'substrate
describe('MSL')
    MSL: 3.2.3 - used components: RealInput, RealOutput
system_info()
     System information
     -OS: Linux
     -Python: 3.10.12
     -Scipy: not installed in the notebook
     -FMPy: 0.3.15
     -FMU by: OpenModelica Compiler OpenModelica 1.21.0
      -FMI: 2.0
      -Type: ME
      -Name: BPL_TEST2_PID.Fedbatch_reg6
     -Generated: 2023-08-22T10:54:51Z
     -MSL: 3.2.3
```

-Description: Bioprocess Library version 2.1.1 -Interaction: FMU-explore for FMPy version 0.9.8

Colab paid products - Cancel contracts here

✓ 0s completed at 20:37

×